组织细胞总胆固醇含量酶法测定试剂盒 Tissue Total Cholesterol (TC) Content Assay Kit 货号: E1015

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1. Ning, H., Jiang, Y., Li, B. et al. CircABCA1 promotes ccRCC by reprogramming cholesterol metabolism and facilitating M2 macrophage polarization through IGF2BP3-mediated stabilization of SCARB1 mRNA. Mol Cancer 24, 199 (2025) (IF:33.9)
2. Liu Y, Dong Y, Jiang Y, et al. Caloric restriction prevents inheritance of polycystic ovary syndrome through oocyte-mediated DNA methylation reprogramming. Cell Metab. 2025 Apr 1;37(4):920-935.e6 (IF:30.9)
3. Wei X, Wu D, Li J, et al. Myeloid beta-arrestin 2 depletion attenuates metabolic dysfunction-associated steatohepatitis via the metabolic reprogramming of macrophages[J]. Cell Metabolism, 2024, 36(10): 2281-2297. e7 (IF:27.7)
4. Xie W, Gan J, Zhou X, et al. Myocardial infarction accelerates the progression of MASH by triggering immunoinflammatory response and induction of periostin[J]. Cell Metabolism, 2024, 36(7): 1634 (IF:27.7)
5. Peng H, Chen B, Wei W, et al. N 6-methyladenosine (m6A) in 18S rRNA promotes fatty acid metabolism and oncogenic transformation[J]. Nature metabolism, 2022, 4(8): 1041-1054 (IF:19.8)
6. Zhang Y, Lin S, Peng J, et al. Amelioration of hepatic steatosis by dietary essential amino acid-induced ubiquitination[J]. Molecular Cell, 2022, 82(8): 1528-1542. e10 (IF:19.3)
7. Fang ZQ, Ruan B, Liu JJ, et al. Notch-triggered maladaptation of liver sinusoidal endothelium aggravates nonalcoholic steatohepatitis through endothelial nitric oxide synthase. Hepatology. 2022 Sep;76(3):742-758

(IF:17.3)

8. Li L, Zhang M, Li J, et al. Cholesterol removal improves performance of a model biomimetic system to co-deliver a photothermal agent and a STING agonist for cancer immunotherapy[J]. Nature communications, 2023, 14(1): 5111 **(IF:16.6)**
9. Liang, M., Ruan, S., Wang, Q. et al. A mitochondria-targeted fluoropolymer nanoparticle with inherent mitophagy-inducing and red fluorescence properties for treatment of atherosclerosis. Nat Commun 16, 9845 (2025) **(IF:15.7)**
10. Guo Y, Zhao M, Bo T, et al. Blocking FSH inhibits hepatic cholesterol biosynthesis and reduces serum cholesterol[J]. Cell research, 2019, 29(2): 151-166 **(IF:15.4)**
11. Wang B, Liang B, Huang Y, et al. Long - Chain Acyl Carnitines Aggravate Polystyrene Nanoplastics - Induced Atherosclerosis by Upregulating MARCO[J]. Advanced Science, 2023, 10(19): 2205876. (2023) **(15.1)**
12. Na S, Fan Y, Chen H L, et al. PPAR α affects hepatic lipid homeostasis by perturbing necroptosis signals in the intestinal epithelium[J]. Acta Pharmaceutica Sinica B, 2024, 14(11): 4858-4873 **(IF:14.7)**
13. Shen, **yue, et al. Hepatic DDAH1 mitigates hepatic steatosis and insulin resistance in obese mice: Involvement of reduced S100A11 expression. Acta Pharmaceutica Sinica B 13.8 (2023): 3352-3364 **(IF:14.5)**
14. Zhang Z, Qin S, Chen Y, et al. Inhibition of NPC1L1 disrupts adaptive responses of drug - tolerant persister cells to chemotherapy[J]. EMBO molecular medicine, 2022, 14(2): e14903 **(IF:14.3)**
15. Lin Y, Ni X, Zhu L, et al. Multi - miRNAs - Mediated Hepatic Lepr Axis Suppression: A Pparg - Dicer1 Pathway - Driven Mechanism in Spermatogenesis for the Intergenerational Transmission of Paternal Metabolic Syndrome[J]. Advanced Science, 2025: 2410831 **(IF:14.1)**
16. Luo X, Huang B, Xu P, et al. The Placenta Regulates Intrauterine Fetal Growth via Exosomal PPAR γ . Adv Sci (Weinh). 2025 Apr;12(15):e2404983 **(IF:14.1)**
17. Xu M, Xu H, Ling YW, et al. Neutrophil extracellular traps-triggered hepatocellular senescence exacerbates lipotoxicity in non-alcoholic steatohepatitis. J Adv Res. 2026 Jan;79:521-534 **(IF:13.0)**
18. Gu X, Hua S, Li C, et al. Konjac glucomannan/alginate polysaccharide hydrogel loaded with phosphatidylcholine nanoparticles for hormone-free cultured fat production from porcine fibroblasts. Carbohydr Polym. 2025 Nov 15;368(Pt 2):124208 **(IF:12.5)**
19. Yin M, Zhang Y, Liu S, et al. The chemokine CCL20 promotes hepatocyte cholesterol deposition during metabolic dysfunction-associated steatohepatitis by regulating OLR1 expression. Metabolism. 2025 Sep;170:156301 **(IF:11.9)**
20. Qiao Y, Zhang Y, Sun C, et al. NSD2 exacerbates metabolic dysfunction-associated steatotic liver disease progression by suppressing TFEB-mediated autophagy-lysosomal pathway. Metabolism. 2026 Jan;174:156416 **(IF:11.9)**
21. Xiao PT, Xie ZS, Kuang YJ, et al. Discovery of a potent FKBP38 agonist that ameliorates HFD-induced hyperlipidemia via mTOR/P70S6K/SREBPs pathway. Acta Pharm Sin B. 2021 Nov;11(11):3542-3552 **(IF:11.9)**
22. Liu R, Li J, Zhang T, et al. Itraconazole suppresses the growth of glioblastoma through induction of autophagy: involvement of abnormal cholesterol trafficking. Autophagy. 2014 Jul;10(7):1241-55 **(IF:11.4)**
23. Liu M, Lu Y, Zhou H, et al. Long-term low-dose exposure to 9-chlorophenanthrene induces liver lipid accumulation via disrupting circadian rhythm. J Hazard Mater. 2025 Sep 5;495:139145 **(IF:11.3)**
24. Feng Y, Xian H, Bai R, et al. Multi-omics insights into 6PPD- and 6PPDQ-induced gut-liver axis disruption and non-alcoholic fatty liver disease progression in zebrafish (Danio rerio). J Hazard Mater. 2025 Sep

5;495:138822 (IF:11.3)

25. Deng CM, Zhang GG, Liu QW, et al. ANO1 Reprograms Cholesterol Metabolism and the Tumor Microenvironment to Promote Cancer Metastasis. *Cancer Res.* 2023 Jun 2;83(11):1851-1865 (IF:11.2)

26. Shen, W., Shao, W., Wang, Q. et al. Dietary diosgenin transcriptionally down-regulated intestinal NPC1L1 expression to prevent cholesterol gallstone formation in mice. *J Biomed Sci* 30, 44 (2023) (IF:11.0)

27. Xin Gu,Shiyuan Hua,Yuqin Huang,et al. κ -Carrageenan/konjac glucomannan composite hydrogel-based 3D porcine cultured meat production[J].*Food Hydrocolloids*, 2024, 151(000) (IF:11.0)

28. Feng J, Ma X, Liu Y, et al.The Role of Human Adiponectin Receptor 1 in 2-Ethylhexyl Diphenyl Phosphate Induced Lipid Metabolic Disruption. *Environ Sci Technol.* 2024 Oct 15;58(41):18190-18201 (IF:10.8)

29. Lin P, Wang Q, Wang Q, et al. Evaluation of the anti-atherosclerotic effect for *Allium macrostemon* Bge. Polysaccharides and structural characterization of its a newly active fructan. *Carbohydr Polym.* 2024 Sep 15;340:122289 (IF:10.7)

30. Chen J, Li H, Yao J, et al.Chronic exposure to PFO4DA and PFO5DoDA, two perfluoroalkyl ether carboxylic acids (PFECAs), suppresses hepatic stress signals and disturbs glucose and lipid metabolism in male mice. *J Hazard Mater.* 2021 Jun 5;411:124963 (IF:10.6)

31. Li Y, Nie JJ, Yang Y, et al. Redox-Unlockable Nanoparticle-Based MST1 Delivery System to Attenuate Hepatic Steatosis via the AMPK/SREBP-1c Signaling Axis. *ACS Appl Mater Interfaces.* 2022 Aug 3;14(30):34328-34341 (IF:10.4)

32. Dai S, Luo M, Jiang T, et al. Dexamethasone as an emerging environmental pollutant: Disruption of cholesterol-dependent synaptogenesis in the hippocampus and subsequent neurobehavioral impacts in offspring. *Environ Int.* 2024 Oct;192:109064 (IF:10.3)

33. Chen, Y., Zhou, J., Tang, C. et al. CD209a regulates metabolic dysfunction-associated steatotic liver disease development through macrophage STAT3 signaling pathway. *Cell Mol Biol Lett* 30, 141 (2025) (IF:10.2)

34. Guo C, Zhao X, Ma R, et al. Silica nanoparticles promoted pro-inflammatory macrophage and foam cell transformation via ROS/PPAR γ /NF- κ B signaling. *Sci Total Environ.* 2023 Jul 10;881:163430 (IF:9.8)

35. Zhang J, Yang W, Zhu Y, et al. Microenvironment-induced programmable nanotherapeutics restore mitochondrial dysfunction for the amelioration of non-alcoholic fatty liver disease. *Acta Biomater.* 2025 Mar 1;194:323-335 (IF:9.6)

36. Du, T., Han, Y., Han, H. et al. The tRNA methyltransferase Mettl1 governs ketogenesis through translational regulation and drives metabolic reprogramming in cardiomyocyte maturation. *Nat Cardiovasc Res* 3, 1438–1453 (2024) (IF:9.4)

37. Fan Y, Sheng S, Guo C, et al.Forkhead Box i2 Transcription Factor Regulates Systemic Energy Metabolism Via Neuropeptide AgRP. *Diabetes.* 2022 Oct 1;71(10):2106-2122 (IF:9.3)

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