活性氧检测试剂盒 (Reactive Oxygen Species Assay Kit) C1300

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1. Liu Xin, Sun Jiao, 2010. Endothelial cells dysfunction induced by silica nanoparticles through oxidative stress via JNK/P53 and NF-κB pathways. *Biomaterials*, 31(32): 8198-8209.
2. Yuan Shaopeng, Fu Yan, Wang Xiaofeng, Shi Hubing, Huang Yujie, Song Xiaomin, Li Ling, Song Nan, Luo Yongzhang, 2008. Voltage-dependent anion channel 1 is involved in endostatin-induced endothelial cell apoptosis. *The FASEB Journal*, 22: 2809-2820.
3. Deng Sanming, Yang Ye, Han Yong, Li Xiaofei, Wang Xiaoping, Li Xueyong, Zhang Zhipei, Wang Yunjie, 2012. UCP2 Inhibits ROS-Mediated Apoptosis in A549 under Hypoxic Conditions. *PLoS ONE*, 7(1): e30714. doi:10.1371/journal.pone.0030714.
4. Zhuo Wei, Song Xiaomin, Zhou Hao, Luo Yongzhang, 2011. Arginine deiminase modulates endothelial tip cells via excessive synthesis of reactive oxygen species. *Biochemical Society Transactions*, 39: 1376-1381.
5. Yue Yang, Zhou Huafeng, Liu Guanlan, Li Yan, Yan Zemin, Duan Mingxing, 2010. The advantages of a novel CoQ10 delivery system in skin photo-protection. *International Journal of Pharmaceutics*, 392(1-2): 57-63.
6. Zhong Z-M, Bai L, Chen J-T, 2009. Advanced Oxidation Protein Products Inhibit Proliferation and Differentiation of Rat Osteoblast-like Cells via NF-κB Pathway. *Cellular Physiology and Biochemistry*, 24(1-2): 105-114.
7. Chen Guozhu, Zhang Xuhui, Zhao Ming, Wang Yan, Cheng Xiang, Wang Di, Xu Yuanji, Du Zhiyan, Yu Xiaodan, 2011. Celastrol targets mitochondrial respiratory chain complex I to induce reactive oxygen species-dependent cytotoxicity in tumor cells. *BMC Cancer*, 11:170. doi:10.1186/1471-2407-11-170.
8. Zhao Hong, Liu GuoLiang, Wang QiuYue, Ding LiYing, Cai Hui, Jiang HaiHong, Xin ZhongQiu, 2007. Effect of ghrelin on human endothelial cells apoptosis induced by high glucose. *Biochemical and Biophysical Research Communications*, 362(3): 677-681.