

蛋白质定量试剂盒 (BCA 法)

P1511

描述: Bicinchoninic acid (BCA)法是近来广为应用的蛋白定量方法。其原理与Lowery法蛋白定量相似，即在碱性环境下蛋白质与Cu²⁺络合并将Cu²⁺还原成Cu¹⁺。BCA与Cu¹⁺结合形成稳定的紫蓝色复合物，在562 nm处有高的光吸收值并与蛋白质浓度成正比，据此可测定蛋白质浓度。与Lowery法相比，BCA蛋白测定方法灵敏度高，操作简单，试剂及其形成的紫蓝色复合物稳定性俱佳，并且受干扰物质影响小。与Bradford法相比，BCA法的显著优点是不受去垢剂的影响。

组成与储存: (1) BCA Reagent 100 ml, 室温保存; (2) Cu Reagent 2 .5ml, 室温保存; (3) BSA standard 4 mg/ml 1 ml, -20℃ 冻存。12 个月有效。可进行 500 次微板(microplate)测定或 100 次 1 ml 比色杯测定。

所需设备: 比色计、酶标仪或微板比色仪，最佳工作波长 562 nm，可在 540-590 nm 之间。

工作溶液(Working Reagent, WR)配制: 将 50 体积 BCA Reagent 与 1 体积 Cu Reagent 混合即为 WR 工作试剂，呈嫩绿色；室温 1 周内稳定。

标准蛋白溶液配制: 用双蒸水、0.9%生理盐水、PBS 或与待测蛋白样品匹配的缓冲液进行倍比稀释: 20 μl 4000 μg/ml BSA + 30 μl 稀释溶液(H₂O/PBS/0.9%NaCl) = 50 μl (BSA=1600 μg/ml)，从中取 25 μl 连续倍比稀释，得到 BSA 标准溶液 1600、800、400、200、100、50、25 μg/ml，各 25 μl。通常，样品蛋白浓度不会太高，也可以预先稀释待测样品，可以省略 1600 μg/ml 标准管而直接从 800 或 1000 μg/ml 开始，能节省标准蛋白用量。

蛋白浓度测定:

蛋白质浓度线性检测范围为 10-2000 μg/ml。标准测定时，用 1 cm 光程玻璃或塑料比色皿，反应终体积 1.1 ml，比色计测定。微板测定时，用 96 孔板，反应终体积 225 μl，用酶标仪、微板比色仪测定。

1. 标准测定时：将 0.05~0.1 ml 标准品或待测样本与 1 ml WR 工作溶液混合。
微板测定时：将 25 μl 标准品或待测样本与 200 μl WR 工作溶液混合。
2. 37°C 反应 30min; 也可 25°C 室温 2 小时或过夜。60°C 30 min 反应可增加检测灵敏度至 5-250μg/ml。
3. 将反应管冷却至室温。测定562 nm (可在540-590 nm之间)光密度(OD)值。
4. 绘制标准曲线。X轴为BSA标准蛋白浓度(mg/ml或μg/ml)，Y轴为各标准管对应的OD562值。用Excel 拟合曲线并计算蛋白浓度。

表 1 标准测定时和微板测定时的加样量和比例

		微板(microplate)测定方案		标准比色杯测定方案	
标号	蛋白浓度 (μg/ml)	标准或待测蛋白体积 (μl)	WR 工作试剂(μl)	标准或待测蛋白体积 (ml)	WR 工作试剂(ml)
1	0	25	200	0.05-0.1	1
2	25	25	200	0.05-0.1	1
3	50	25	200	0.05-0.1	1
4	100	25	200	0.05-0.1	1
5	200	25	200	0.05-0.1	1
6	400	25	200	0.05-0.1	1
7	800	25	200	0.05-0.1	1
8	1600	25	200	0.05-0.1	1
待测样品		25	200	0.05-0.1	1

注意事项

1. 37°C 30min 或 25°C 室温反应 2 小时对测量较为便利，但严格来讲此时反应尚未达到终点，通常每 10 min OD562 值升高约 2.3%。然而，通常在 10min 内可以测定 30 管并不会明显影响测定精度。
2. BCA 法检测范围为 20-2000 μg/ml。检测 0.5~10 μg/ml 高度稀释蛋白样品应采用 Bradford 法蛋白质

定量试剂盒(# P1510)。60°C 30 min 反应可增加检测灵敏度至 5~250μg/ml。

3. BCA法检测样品中含有脂类物质时光吸收值会偏高。样品中EDTA或葡萄糖浓度大于10 mM不能使用BCA方法，葡萄糖浓度大于10 mM时可用改良Lowry法蛋白定量试剂盒 #P1512, EDTA大于10 mM的样品可用Bradford法蛋白质定量试剂盒(# P1510)。另外，蛋白质样品经液体样品蛋白抽提试剂(#P1255)沉淀后，可彻底去除干扰BCA法、Bradford法、和Lowry法蛋白测定的物质。
4. 欲使测量能耐受下面表 2 所提示的最大干扰物质浓度，并保持测量精度，应在蛋白标准管中加入相应浓度的干扰物质，但会给操作带来不便。
5. 可测量吸附于固相支持物乳酶标板、琼脂糖、亲和层系凝胶上的蛋白。
6. 每次测定应该重新测定并制作标准曲线。

参考文献:

Smith P et al, 1995, Measurment of protein using bicinchiconic acid, Anal. Biochem. 150, 76-85

表2 BCA法物质干扰及耐受的最大浓度

Buffer Systems	Sodium phosphate 25 mM
Bicine, pH 8.4 20 mM	Sucrose 40%
Bis-Tris, pH 6.5 33 mM	Sodium ortho-Vanadate in PBS, pH 7.2, 1 mM
Calcium chloride in TBS, pH 7.2 10 mM	Urea 3 M
CHES, pH 9.0 100 mM	Chelating agents
Cobalt chloride in TBS, pH 7.2 0.8 M	EDTA 10 mM
Ferric chloride in TBS, pH 7.2 10 mM	EGTA,any level, not compatible
HEPES 100 mM	Sodium citrate 200 mM
MOPS, pH 7.2 100 mM	Detergents
Nickel chloride in TBS 10 mM	Brij-35 5%
PBS; no interference	Brij-52 1%
NaCl (0.15 M), pH 7.2, no interference	CHAPS 5%
PIPES, pH 6.8 100 mM	CHAPSO 5%
Sodium acetate, pH 4.8 200 mM	Deoxycholic acid 5%
Sodium citrate, pH 4.8 or pH 6.4 200 mM	Nonidet P-40 (Igepal CA-630) 5%
Tricine, pH 8.0 25 mM	SDS 5%
Triethanolamine, pH 7.8 25 mM	Span 20 1%
Tris 250 mM	Triton X-100 5%
TBS buffer, no interference	Triton X-114 1%
1 x SDS-PAGE loading buffer, no interference	Tween-20 5%
Zinc chloride (10 mM) in TBS, pH 7.2, 10 mM	Tween-60 5%
Buffer Additives	Tween-80 5%
Ammonium sulfate 1.5 mM	Zwittergents 1%
Aprotinin 10 mg/L	Reducing & Thiol Containing Agents
Glucose 10 mM	Dithioerythritol (DTE) 1 mM
Glycerol 10%	Dithiothreitol (DTT) 1 mM
Guanidine•HCl 4 M	2-Mercaptoethanol 1 mM
HCl 100 mM	Tributyl Phosphine 0.01%
Imidazole 50 mM	Solvents
Leupeptin 10 mg/L	Acetone 10%
PMSF 1 mM	Acetonitrile 10%
Sodium azide 0.20%	DMF 10%
Sodium bicarbonate 100 mM	DMSO 10%
Sodium chloride 1 M	Ethanol 10%
Sodium hydroxide 100 mM	Methanol 10%

目前，使用普利莱 BCA 蛋白定量试剂盒发表的 SCI 文章已超过百篇，部分文章如下，供参考。

- 1、Gong X, Chen Z, Liu Y, et al. Gene expression profiling of the paracrine effects of uterine natural killer cells on human endometrial epithelial cells[J]. International journal of endocrinology, 2014, 2014.
- 2、Wang W, Zhao J, Yang G, et al. Heat shock protein 47 is up-regulated in conjunctival bleb scarring after filtration surgery in rats[J]. Graefe's Archive for Clinical and Experimental Ophthalmology, 2015, 253(1): 91-97.
- 3、Hou G, Li C, Liu S, et al. Long-term treatment with EXf, a peptide analog of Exendin-4, improves β -cell function and survival in diabetic KKAY mice[J]. Peptides, 2013, 40: 123-132.
- 4、Zhang Q H, Li J C, Dong N, et al. Burn injury induces gelsolin expression and cleavage in the brain of mice[J]. Neuroscience, 2013, 228: 60-72.
- 5、Zhang Y, Liu X, Zuo T, et al. Tetramethylpyrazine reverses multidrug resistance in breast cancer cells through regulating the expression and function of P-glycoprotein[J]. Medical Oncology, 2012, 29(2): 534-538.
- 6、Mo Y, Yao H, Song H, et al. Alteration of Behavioral Changes and Hippocampus Galanin Expression in Chronic Unpredictable Mild Stress-Induced Depression Rats and Effect of Electroacupuncture Treatment[J]. Evidence-Based Complementary and Alternative Medicine, 2014, 2014.
- 7、Liu C, Huang H, Wang C, et al. Involvement of ephrin receptor A4 in pancreatic cancer cell motility and invasion[J]. Oncology letters, 2014, 7(6): 2165-2169.
- 8、Wei N, Liu G T, Chen X G, et al. H1, a derivative of Tetrandrine, exerts anti-MDR activity by initiating intrinsic apoptosis pathway and inhibiting the activation of Erk1/2 and Akt1/2[J]. Biochemical pharmacology, 2011, 82(11): 1593-1603.
- 9、Yang Y, Li C, Xiang X, et al. Ursolic acid prevents endoplasmic reticulum stress-mediated apoptosis induced by heat stress in mouse cardiac myocytes[J]. Journal of molecular and cellular cardiology, 2014, 67: 103-111.
- 10、Lin F S, Shen S Q, Chen Z B, et al. 17 β -Estradiol attenuates reduced-size hepatic ischemia/reperfusion injury by inhibition apoptosis via mitochondrial pathway in rats[J]. Shock, 2012, 37(2): 183-190.
- 11、Li H, Chen Q, Liu F, et al. Unfolded protein response and activated degradative pathways regulation in GNE myopathy[J]. PloS one, 2013, 8(3): e58116.
- 12、Ma H, You G P, Zhang X P, et al. A novel role of globular adiponectin in treatment with HFD/STZ induced T2DM combined with NAFLD rats[J]. The Scientific World Journal, 2014, 2014.
- 13、Li X, Liu J, Lin L, et al. Traditional Chinese Medicine Shuang Shen Ning Xin Attenuates Myocardial Ischemia/Reperfusion Injury by Preserving of Mitochondrial Function[J]. Evidence-Based Complementary and Alternative Medicine, 2014, 2014.
- 14、Ge Z H, Wang Z X, Yu T L, et al. Morphine improved the antitumor effects on MCF-7 cells in combination with 5-Fluorouracil[J]. Biomedicine & Pharmacotherapy, 2014, 68(3): 299-305.
- 15、Lin C, Wang Y, Lai Y, et al. Incorporation of carboxylation multiwalled carbon nanotubes into biodegradable poly (lactic-co-glycolic acid) for bone tissue engineering[J]. Colloids and Surfaces B: Biointerfaces, 2011, 83(2): 367-375.
- 16、Liang Y, Li Z, Mo N, et al. Isoflurane Preconditioning Ameliorates Renal Ischemia-Reperfusion Injury through Antiinflammatory and Antiapoptotic Actions in Rats[J]. Biological & pharmaceutical bulletin, 2014, 37(10): 1599-1605.
- 17、Zhao X, Gao M, He J, et al. Perilipin1 Deficiency in Whole Body or Bone Marrow-Derived Cells Attenuates Lesions in Atherosclerosis-Prone Mice[J]. PloS one, 2015, 10(4).
- 18、He B, Lin P, Jia Z, et al. The transport mechanisms of polymer nanoparticles in Caco-2 epithelial cells[J]. Biomaterials, 2013, 34(25): 6082-6098.
- 19、He B, Lin P, Jia Z, et al. The transport mechanisms of polymer nanoparticles in Caco-2 epithelial cells[J]. Biomaterials, 2013, 34(25): 6082-6098.
- 20、Zhang M, Xu J J, Zhou R L, et al. cAMP responsive element binding protein-1 is a transcription factor of lysosomal-associated protein transmembrane-4 beta in human breast cancer cells[J]. PloS one, 2013, 8(2): e57520.
- 21、Dang Y M, Huang G, Chen Y R, et al. Sulforaphane inhibits the proliferation of the BIU87 bladder cancer cell line via IGFBP-3 elevation[J]. Asian Pac J Cancer Prev, 2014, 15: 1517-20.